

Inactivity Monitor

Working with the IntelliLink Executive, the PIR can look for movement within a defined time. If the user has not moved within that period, due to a fall or feeling unwell, the carephone will send an alarm to the carer/monitoring centre.

Intruder Monitor

If the user leaves their property to stay with friends or visit a day centre, pressing a button on the carephone and the portable button together will set the system into 'away' mode. The 'away' mode will detect intruders, giving peace of mind to the user.



Technical Data	
Model Number	FG159048
Warranty	1 years warranty
Power Requirements	Alkaline or Lithium Batteries. Typical battery life 2 years depending on usage.
Size	Height 110mm x Width 65mm x Depth 75mm

Chubb Community Care Shadsworth Road, Blackburn
BB1 2PR Tel: 01254 688 774 Fax: 01254 696 460

